

TBC1D29-TBC1 domain family, member 29 Gene

Reference: PTXBC096718

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TBC1 domain family member 29		
For research use only.			

GENE DESCRIPTION	TBC1 domain family, member 29
GENE ID	26083
ORIGIN SPECIES	Human
ACCESSIONS	BC096718
VECTOR	pDONR223
NCBI SYMBOL	TBC1D29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcatctggacaaggaaggtctatgcacacaggggtctcattcagctggcttccgggtgctgaatgatggatctctcgggctcaccctgcctgtgggatatgtatttgcggagggaacagatgtgatgctgataacaagcattgccttaagggtcaaaggagcctcta
tgaagaaactaacaaggaacatgggacctgccacccagggccctcaaggcactggaagagccagaccatttgtgagagcctccactctccctgcaagcgtgacagcctcagagagcagagagggccctcactcctgcagactctccaagggtgccaggacaacaag
ccttgagccggggagacaagggaatcagtggtctcattatctctaccatctctaccatctcggagaggagatgtggcagataatagggttaa