

TP53TG3-TP53 target 3 Gene

Reference: PTXBC096733

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P53TG3; TP53TG3A; TP53TG3E; TP53TG3F; TP53-target gene 3 protein; TP53-inducible gene 3 protein; p53 target gene 3; TP53 target 3		

For research use only.

GENE DESCRIPTION	TP53 target 3
GENE ID	24150
ORIGIN SPECIES	Human
ACCESSIONS	BC096733
VECTOR	pDONR223
NCBI SYMBOL	TP53TG3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgcgcctcaccctgcacgtccagcccgagccagctggcctcctagaccctctgccctgcgaccaacagccgggagcggaccagacaccagaactccggacgggtgaagacgggttcgctccctgtccgccttcgcagcccagcagtttgcctcgcggagaggagccttgctgtttccaaatctcctcgtgaagagacattggagctagggcggcgtagtttcacctggctgtcacaatctctgcacccagcagcggcgtggaccccgggcctcgtcgtgtgtccatctag