

PART1-prostate androgen-regulated transcript 1 Gene

Reference: PTXBC098111

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCRNA00206; prostate androgen-regulated transcript 1 (non-protein coding)		
For research use only.			

GENE DESCRIPTION	prostate androgen-regulated transcript 1
GENE ID	25859
ORIGIN SPECIES	Human
ACCESSIONS	BC098111
VECTOR	pDONR223
NCBI SYMBOL	PART1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaatgacgctatttagacagaacatccaaggccgtgtcagaactcaattacgactacatatgcattaaggcaggaactggcaggcctcagggtacgccaactataggactcgtgcttctcgtacgctggcctataatctatgaaactgagctccagagccagccaatcacttag