

ABCA11P-ATP-binding cassette, sub-family A (ABC1), member 11 (pseudogene) Gene

Reference: PTXBC098116

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ATP-binding cassette, sub-family A (ABC1), member 11 (pseudogene)
GENE ID	79963
ORIGIN SPECIES	Human
ACCESSIONS	BC098116
VECTOR	pDONR223
NCBI SYMBOL	ABCA11P
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccctgaggggcagcagcaaatgtggcagatactcaggctaccattaaaaaccaggagagggggccctcttgaccaccattacatgcagaggttaagtctctgtgtgaccgtgtggccatcatggtgtcaggacgtaagggtattggtccattcaacagctgaaaagtgtg
gtaagattattactag