

ULBP3-UL16 binding protein 3 Gene

Reference: PTXBC098127

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	N2DL-3; NKG2DL3; RAET1N; NKG2D ligand 3; ALCAN-gamma; retinoic acid early transcript 1N; UL16 binding protein 3		

For research use only.

GENE DESCRIPTION	UL16 binding protein 3
GENE ID	79465
ORIGIN SPECIES	Human
ACCESSIONS	BC098127
VECTOR	pDONR223
NCBI SYMBOL	ULBP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggcagcggcgcagccccgcgatcctccgcgcctcgcgattctccgtacctgctattcgaactgtccggacggggcggccgacgctcactctctgtataactcaccatcattcattgccagacatgggcaacagtgtgtgaggtccagaccaggtgatcagaaga
attttctcctatgactgtgagtgacaaggtcttatctatgggtcacctagaagagcagctgtatgccacagatgctgggaaaacaactggaaatgctgagagagtgggcgagagctcagactggaactggctgacactgagctggaggattcacaccagtgaccctcac
gctgcaggtcaggatgctgtgagtgtaagccgatggataatccgtggatctggcagttcagctcgtatggacggaagtctcctctcttctgactcaaacacagaaaatggacagtggttcacgctggagccagggcggatgaaagagaagtgaggagaaggtatgactgaccac
ctcttcaagatggtcctcaatgagagactgcaagagctggcttagggacttctgatgcacaggaagaagaggtggaaccacagaccaccaccatgccccaggttagctcaacccaagccatagccaccacctcagtcctggagttctctatcctctgcttcatctcc
ctggcatctga
```