



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

VCX3A-variable charge, X-linked 3A Gene

Reference: PTXBC098149

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VCX-8r; VCX-A; VCX3; VCX8R; VCXA; variable charge X-linked protein 3; variable charge protein on X with eight repeats; variably charged X-A; variably charged protein X-A; variable charge, X-linked 3A		

For research use only.

GENE DESCRIPTION	variable charge, X-linked 3A
GENE ID	51481
ORIGIN SPECIES	Human
ACCESSIONS	BC098149
VECTOR	pDONR223
NCBI SYMBOL	VCX3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttccaaagccgagagcgtcgggaccctcggccaaagccacggaggcaggaagaagggaagtcctctctcagccgagcccgatgaccccgaagaagaagactaccaagctggcccaagaagggaagaacagttctgtagagggaagacccgggaagaagaagggtctgcgaca
aagatggccggccgtgacggacatcctggggcggagagccggccgagcggccacccggccgagcagcagccagccagcagcagtgagtcctctcagcagcagctctccgcggcgaagagcagatgagcagcagggagccagccagcagccctgagctcagggagagcgtgctggaggaaccac
tctgacgagagccgaggtgctggaagaaccactgagtcagcagagccagctgaggaaccactcgtgacgagagcaggtgctcagcaggaaccactgagtcagcagagccagcgtgagggaagaaccactgagtcagcagagcagagatggaagaactcagatgatg

More info: Please contact us for any question: contact@proteogenix.fr