

KCNQ1DN-KCNQ1 downstream neighbor Gene

Reference: PTXBC098159

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BWRT; HSA404617; KCNQ1 downstream neighbor (non-protein coding)		

For research use only.

GENE DESCRIPTION	KCNQ1 downstream neighbor
GENE ID	55539
ORIGIN SPECIES	Human
ACCESSIONS	BC098159
VECTOR	pDONR223
NCBI SYMBOL	KCNQ1DN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcaggaagtgtcaggaccactgccgagaccaactccccatgccaccaccagggtgcgcctggactcctggaagggtgcgctcagggcagccatcaaaggctccaagagcaagaggcaaggagagaagtgtcctactttaacggccagcctcagtggtcggc
cctgtcactctgccacctcagagagatga