

PROK2-prokineticin 2 Gene

Reference: PTXBC098162

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BV8; HH4; KAL4; MIT1; PK2; prokineticin-2; protein Bv8 homolog; prokineticin 2		
For research use only.			

GENE DESCRIPTION	prokineticin 2
GENE ID	60675
ORIGIN SPECIES	Human
ACCESSIONS	BC098162
VECTOR	pDONR223
NCBI SYMBOL	PROK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggagcctgtgtgcgccccactctgtctcttctgtgtgcgcgcgtgtgtcagccccgcgtggggacgcgcgcgtgtatcaccgggcttgtgacaaggactccaatgtgtggaggcatgtgtgtgtgtcagtcctgtgggtcaagagcataaggattgcacacctat
gggcaaacctgggagacagctgccatccactgactcgtaaaaacaattttggaatggaaggcaggaaagaagaaagggaagagaagcaaaaggaaaaaggaggtccatttttggcgaggatgcatcacactgcccatgtctgccaggcttggcctgtttacggactcatttaac
cgatttattgtttagcccaaaagtaa