

LOC729355-similar to TP53TG3 protein Gene

Reference: PTXBC098165

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TP53TG3E; TP53TG3F; TP53-target gene 3 protein; TP53-inducible gene 3 protein; TP53 target 3B		

For research use only.

GENE DESCRIPTION	similar to TP53TG3 protein
GENE ID	729355
ORIGIN SPECIES	Human
ACCESSIONS	BC098165
VECTOR	pDONR223
NCBI SYMBOL	LOC729355
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgcgcctcaccctgcacccccagccagcagctggcctcctagaccctctgccctgcgaccaacagccgggagcggaccagaccagaactccggacgggtgaagacgggtccgcctctcccgcccttcgcagccagcagtttcgccctgcggagaggagcctt
gctgtttccaaatctcctcctgctgaagagacattggagctagggcggcgtagtttcacctggtaattgtgacaccctgtctcctcgagctgcaggcttttatgcttgcattgttcgaagtttgataccttgcagatcaacaaggccgggtggcctctcactgcctcgcggcagggtgtcaagtttt
tcagggttaa