

LENEP-lens epithelial protein Gene

Reference: PTXBC098246

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	lens epithelial cell protein LEP503; lens epithelial protein		

For research use only.

GENE DESCRIPTION	lens epithelial protein
GENE ID	55891
ORIGIN SPECIES	Human
ACCESSIONS	BC098246
VECTOR	pDONR223
NCBI SYMBOL	LENEP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagccccggacacagccccctagcccaaacctacccttctctcggaggggccctcgagacactgggctcgggtgcctgtcattaagatgggcacaggggtggagggtccagcggaccctgaagggaagtcgcctacatcctctctgctgctgggtgatcaagggaactgctggattaa