



Reference: PTXBC098248

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Unc-93A; dJ366N23.1; dJ366N23.2; protein unc-93 homolog A; dJ366N23.1 (putative C. elegans UNC-93 (protein 1, C46F11.1) LIKE); dJ366N23.2 (putative C. elegans UNC-93 (protein 1, C46F11.1) C-terminal LIKE); unc93 homolog A; unc-93 homolog A (C. elegans)		

For research use only.

GENE DESCRIPTION	unc-93 homolog A (<i>C. elegans</i>)
GENE ID	54346
ORIGIN SPECIES	Human
ACCESSIONS	BC098248
VECTOR	pDONR223
NCBI SYMBOL	UNC93A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]