

NKAIN1-Na⁺/K⁺ transporting ATPase interacting 1 Gene

Reference: PTXBC098249

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM77C; sodium/potassium-transporting ATPase subunit beta-1-interacting protein 1; Na(+)/K(+)-transporting ATPase subunit beta-1-interacting protein 1; Na ⁺ /K ⁺ transporting ATPase interacting 1; family with sequence similarity 77, member C; Sodium/potassiu		

For research use only.

GENE DESCRIPTION	Na ⁺ /K ⁺ transporting ATPase interacting 1
GENE ID	79570
ORIGIN SPECIES	Human
ACCESSIONS	BC098249
VECTOR	pDONR223
NCBI SYMBOL	NKAIN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagtcacctgggcatcttggcaccgtgcagtagccggtacccggtacccctgtatgcagcctggctgctggtgggtggaatgcattatcatctgcttacttgagggtggacagctgtccaggaccgggacttcatcatgacctcaacacatccctgcaccgctcc
tgggtgatggagaatggccaggctgacctggtgacacctgtctgaactcccgctggtctggaggaccaccatgtcatctctgtcactggctgctgctgactacccctacattgaagccctcagcagcgccctgcagatcttctggcactgttcggcttctgttcgctgtctacgtga
gcaaagtgttctggaggaggaggacagcttgcattcatggcggcttgcactcctacggataccaggcgccccagaagacgtgcattacagctgcagcctctgtacacgtcgggtag