

NANOGP8-Nanog homeobox pseudogene 8 Gene

Reference: PTXBC098275

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	Nanog homeobox pseudogene 8
GENE ID	388112
ORIGIN SPECIES	Human
ACCESSIONS	BC098275
VECTOR	pDONR223
NCBI SYMBOL	NANOGP8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtggtgacccagctgtcccaaaagctgccttgccttgaagcatccgactgtaagaatcttcacatgacgtgatttgggctggaagaaactatccatccttgcaaatgtcttctgctgagatgccctcacacagagactgtctctctctctccatggatctgcttattcagga
cagccctgattcttccaccagtcctccaaaggcaaacacccactcttcagagaaatagtgctgcaaaaaaggagacaaaggtcccggtcaagaaacagaagaccagaactgtgtctcttccaccagctgtgtgtactcaatgatagattcagagacagaaatccctcagctccagcag
atgcaagaactctcacaatcctgaacctcagctacaacagggtgaagacctggtccagaaccagagaatgaaatctaagaggtggcagaaaaacaactggccgaagaatagcaatggtgtgacgcagaaggcctcagcactacccagcctctactctctaccaccagg
atgcctgtgaaaccgactgggaaccttcaatgtgggagcaaccagacctggaacaattcaacctggagcaaccagaccagaacatccagtcctggagcaaccactctggaaacactcagacctggtgacccaatcctggaaacatcaggcctggaacagtccttataaactgtgg
agaggaatctctgcagtcctgcatgcacttcagccaaattctcctgcagtgactggaggctgccttggaaagctgctggggaaggccttaattgtaatacagcagaccactaggtattttagtactccacaaccatggatttattcctaaactactccatgaacatgcaacctgaagacgtgtg
a