



Reference: PTXBC098278

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CSD; PCAP; P-selectin cytoplasmic tail-associated protein; cysteine sulfinic acid decarboxylase-related protein; cysteine-sulfinatase; sulfinoalanine decarboxylase		

For research use only.

GENE DESCRIPTION	cysteine sulfinic acid decarboxylase
GENE ID	51380
ORIGIN SPECIES	Human
ACCESSIONS	BC098278
VECTOR	pDONR223
NCBI SYMBOL	CSAD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]