



Reference: PTXBC098298

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tubulin monoglycyclase TTLL3; HOTTL; tubulin tyrosine ligase-like family, member 3; tubulin--tyrosine ligase-like protein 3; tubulin tyrosine ligase like 3		

For research use only.

GENE DESCRIPTION	tubulin tyrosine ligase-like family, member 3
GENE ID	26140
ORIGIN SPECIES	Human
ACCESSIONS	BC098298
VECTOR	pDONR223
NCBI SYMBOL	TTLL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr