

IL1F9-interleukin 1 family, member 9 Gene

Reference: PTXBC098337

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL1F9; IL-1F9; IL-1H1; IL-1RP2; IL1E; IL1H1; IL1RP2; interleukin-36 gamma; IL-1 related protein 2; IL-1(EPSILON); IL-1-epsilon; interleukin 1-related protein 2; interleukin-1 epsilon; interleukin-1 family member 9; interleukin-1 homolog 1; interleukin 36,		

For research use only.

GENE DESCRIPTION	interleukin 1 family, member 9
GENE ID	56300
ORIGIN SPECIES	Human
ACCESSIONS	BC098337
VECTOR	pDONR223
NCBI SYMBOL	IL1F9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaggcactccaggagacgtgatggaggaggccgtctatcaatcaatgtgtaaacctattactgggactattaatgattgaatcagcaagtgtggaccctcagggtcagaaccttggcaggtccacgaagtgcacagtgacagtgtagcccgagtcactgtgctgtatcacatgca
agtatccagaggctcttgagcaaggcagaggatccatttatgtggaatccagaatccagaaatgtgtgtatgtgagaaggttgagaacagccacattgcagctaaaagagcagaagatcatggaatctgtatggccaacccgagcccgtagaaccttcctttctaccgtgcca
agactggtaggacctccaccctgagtcgtggcctcccgactgggttcattgcctcctccaagagagaccagcccatctctgacttcagaacttgggaagtcatacaacactgccttgaattaaataataatgactga