

IFRG15-interferon responsive gene 15 Gene

Reference: PTXBC098348

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFRG15; LULL1; NET9; torsin-1A-interacting protein 2; interferon alpha responsive protein; 15 kDa interferon-responsive protein; interferon responsive gene 15; lumenal domain like LAP1; torsin A interacting protein 2; torsin 1A interacting protein 2		

For research use only.

GENE DESCRIPTION	interferon responsive gene 15
GENE ID	64163
ORIGIN SPECIES	Human
ACCESSIONS	BC098348
VECTOR	pDONR223
NCBI SYMBOL	IFRG15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttcagataattcacattgccctgattgtggacaacagtgggtccctagtttagaactaggccactgggtgtaccaaactgaactgttgaataatgaattaccaggtattcttagaccgtattacacagactgattattgtcctgagtggttatcctgataatcctgctaataagaagcctgtgtcttc
cttggtctttccacttgagtggtgctcccgagaatctcaccagatggacctttgagaagcttgcacatcatttctctgggtcctccactggttagaaaaagaatacatgactctcgagtagctgggttttaaccctgcattacagttaatcttgaccagaacagataaaaccttaacaaaaaactg
ggcctaaacaaatag