



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MC3R-melanocortin 3 receptor Gene

Reference: PTXBC098351

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

BMIQ9; MC3; MC3-R; OB20; OQTL; melanocortin receptor 3; obesity quantitative trait locus; melanocortin 3 receptor

For research use only.

GENE DESCRIPTION	melanocortin 3 receptor
GENE ID	4159
ORIGIN SPECIES	Human
ACCESSIONS	BC098351
VECTOR	pDONR223
NCBI SYMBOL	MC3R
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcatccaaaagcagctatctgaggggagatttcttctctctgagcagcagcagcttctacggaccctgctggagcccccagctcgatcagccctctctgacagcaatgaatgctctgctgctcctcctctgttcagcaaacactgcctaatggctcggagcaacctccaaagcccttcttcagcaaacagcagcagcagcctctatgtagcaggctctacacacggccagcgtttctctctctgggcatctgcagctgctggaaaacatctcgtttatctgcgcgtgctcaggaagcaacctcactcccgatgactcttctctctgacacccctggcgctgcagcatgctgaagtgtgtccaatgcccggagaccatcatgatcgcagctctccagcagctactcagcctcaggagcagctttatcagcagctggacaacatctcagatcctcagatgactctctctctggctcctcatcagcactctcgttccctcagctcgcagcaggtacgtcaccatcttttgcgcgtccgctcaccacagcatcatgaccgtgaggaaagccctcaccctgcatctgtggccatctgggtctgctgcgcggctctctggcgtgttcatcgtctacgtcggagcaaaaatggctatctgtgtgctctcatcaccatgttcttcccatgatctctcatgcccaccctctacgtgcacatgtttctcttggcgcgtcagctcagctcagcagcatagcagcctgaccactcgcagcggggggcagctaccatacactcttctcctggcgcttctcatctctctgctggcggccctcttccctcactggtctctcatcactcactcactcagctgcccaaccaacctcatcagctctacacctgcgccacttaccagctctcatcattcagcctcgtctacgacctcactcgttctcggagcgtggaattgctggctgacacagcagctctgggtagat

More info: Please contact us for any question: contact@proteogenix.fr