



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

IFNA6-interferon, alpha 6 Gene

Reference: PTXBC098357

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFN-alphaK; interferon alpha-6; IFN-alpha-6; interferon alpha-54; interferon alpha-K; IeIF K; interferon alpha 6		

For research use only.

GENE DESCRIPTION	interferon, alpha 6
GENE ID	3443
ORIGIN SPECIES	Human
ACCESSIONS	BC098357
VECTOR	pDONR223
NCBI SYMBOL	IFNA6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttcgcttcttcttactgatggccctgggtggctgactcagaatcaagctctctctggactgtgacttgcctcagaccacagcctgggttcacagaggaccatgatgctctggcacaatgaggagaatctctcttctctctgtctgaaggacagacatgacttcagatttcccc
aggaggaggtttgatgccaacaggctccagaaggctgaagccatctctgtctccatgatgggtgaltcagacagacttcaacctcttcagcacaaggagactcatctgttctgttggatgagagggctcttagacaacctatatactgaactttaccacagctgaatgactctggaagcctgtgtgatg
caggagctgtgtgggtggaggagactccccctgatgaatgaggactctatctgctgtgtgagaaaatacttccaaagaatcatctctacctgcagagagaaaaagtacagccctgtgtcctgggaggtgttcagagcagaatatcatgatctctcttcataagaacctgcagaagaagtta
aggaggagaagataa

More info: Please contact us for any question: contact@proteogenix.fr