

FAM106A-family with sequence similarity 106, member A Gene

Reference: PTXBC098363

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	family with sequence similarity 106, member A		
For research use only.			

GENE DESCRIPTION	family with sequence similarity 106, member A
GENE ID	80039
ORIGIN SPECIES	Human
ACCESSIONS	BC098363
VECTOR	pDONR223
NCBI SYMBOL	FAM106A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtaccgtctactatgttcctcgtgcaccttcactcagcacaacaggctgcactgcctccgaacacatcaettgaactcgtcctctgctcctttgtccatctgaatcacctctccacattctgacctgtaactcctgatcctcactcatgaggcagttcgtctcttttcatttagctttcccc
ggggtacactgtcaatagctactgtctgatgtcatcagtcagcacttcacagaggcaattatgtccacagagcttttgcttaattactgtcactcctctttacatgtgtgcatatgcatactcttccctaagagactggaaatcatgattcctctggagcagtggtttccataagtgaccaa
cccacagaccagtgcaaaactggctgcaaaagagcttcctctgaggaactattagaatgcagattcttgcactgcattgggagaggagatctataaactctgggtgtaataggcacagaacgttag