

CCL16-chemokine (C-C motif) ligand 16 Gene

Reference: PTXBC099662

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CKb12; HCC-4; ILINCK; LCC-1; LEC; LMC; Mtn-1; NCC-4; NCC4; SCYA16; SCYL4; C-C motif chemokine 16; IL-10-inducible chemokine; chemokine (C-C motif) ligand 16; chemokine LEC; liver CC chemokine-1; liver-expressed chemokine; lymphocyte and monocyte chemoattr		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 16
GENE ID	6360
ORIGIN SPECIES	Human
ACCESSIONS	BC099662
VECTOR	pDONR223
NCBI SYMBOL	CCL16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggtctccgaggtgacctgtctctcttgcctcatcttattcattctcggtctcgcagccagccaaaagttcctgagtggtgaacaccccatccacctgctgctgaagtattatgagaaagtgttccaaggagactagtgtgtggatagacagaaaggccctcaactgtcacctgccagcaatcatcttcgtacccaagagggaaccgagaagtctgcaccaacccaatgacgactgggtccaagagtacatcaaggatccaacctactttgtgctaccaggaactgtccacgggttaaaattattacagcaagaatggtcaaccccgactctcaactccagtgga