



Reference: PTXBC099703

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BTN5; PRO2801; erythroid membrane-associated protein; Radin blood group (Rd); Scianna blood group (Sc); radin blood group antigen; scianna blood group antigen; erythroblast membrane associated protein (Scianna blood group)		

For research use only.

GENE DESCRIPTION	erythroblast membrane-associated protein (Scianna blood group)
GENE ID	114625
ORIGIN SPECIES	Human
ACCESSIONS	BC099703
VECTOR	pDONR223
NCBI SYMBOL	ERMAP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]