

C12orf48-chromosome 12 open reading frame 48 Gene

Reference: PTXBC099734

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf48; AROM; PARI; PCNA-interacting partner; PARP1 binding protein		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 48
GENE ID	55010
ORIGIN SPECIES	Human
ACCESSIONS	BC099734
VECTOR	pDONR223
NCBI SYMBOL	C12orf48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgtgactgaccattatgaggacgttaggaagatttatgatatttctgaagaacagtaataatgttagatctgattgatgtttatcaaaatgtagggctttgacttctaatttgaaaattataacagatctctagtcactactggattttctgtctggcaaacagtagcagtagtgat
gaaactgatctttctataccaacatcaccaacaagtaataacaaccgtgataatgaaaagggtgcagctgctagcaaggaaaattatctttcatatttaaatctgtagtgaaatcaaagaatgacctggctgtgcttatattctcaatatctctgataaggactaggagaagagccttcactgat
ttgaaacatgctgctcgagagaacaatgtctatcttttgggtggccacgtcttttattagaacaatagagcttggagggaaggatagcaccaccaccatcagatcctttaaggacacatgtaagggtgattgtctaattttatttaatttcattgacaaattagatgagattcttggagaatacaca
accgaagggtgttaaatccatctgttgaagatcaacaattggaacgagttttgaaatgttcacatctggacagaagtaa