



Reference: PTXBC099916

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RBM35A; RMB35A; epithelial splicing regulatory protein 1; RNA-binding motif protein 35A; RNA-binding protein 35A		

For research use only.

GENE DESCRIPTION	RNA binding motif protein 35A
GENE ID	54845
ORIGIN SPECIES	Human
ACCESSIONS	BC099916
VECTOR	pDONR223
NCBI SYMBOL	RBM35A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]