

ZNF663-zinc finger protein 663 Gene

Reference: PTXBC100775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	zinc finger protein 663
GENE ID	100130934
ORIGIN SPECIES	Human
ACCESSIONS	BC100775
VECTOR	pDONR223
NCBI SYMBOL	ZNF663
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtacacaggagaaaaacctgatgagtgtaaaagaaatgagaaagccttgaagaagttataccatactcaaaatgagagaactcatgcaggagagaaagtctatgattgtaagaaatgtgggaaattcttcatgaaaaagcatgccttactcaacatcagagaaccacacacagggaa
accagcaagtgtaatgatagtgaaagagtcctgaaggaaatcttcctcactcccaatcagagaattaaaacaagattaaaacaagagagaaaaactgtaaggtataaatgtgagaaaccttcttggagaagtga