



Reference: PTXBC100809

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GLUT10; GLUT11; solute carrier family 2, facilitated glucose transporter member 11; facilitative glucose transporter GLUT11; glucose transporter protein 10; glucose transporter protein 11; glucose transporter type 10; glucose transporter type 11; glucose		

For research use only.

GENE DESCRIPTION	solute carrier family 2 (facilitated glucose transporter), member 11
GENE ID	66035
ORIGIN SPECIES	Human
ACCESSIONS	BC100809
VECTOR	pDONR223
NCBI SYMBOL	SLC2A11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]