



Reference: PTXBC100826

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	synaptic vesicle glycoprotein 2C; synaptic vesicle protein 2C		
For research use only.			

GENE DESCRIPTION	synaptic vesicle glycoprotein 2C
GENE ID	22987
ORIGIN SPECIES	Human
ACCESSIONS	BC100826
VECTOR	pDONR223
NCBI SYMBOL	SV2C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr