

DEFB106A-defensin, beta 106A Gene

Reference: PTXBC100844

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BD-6; DEFB-6; DEFB106; beta-defensin 106; beta-defensin 6; defensin, beta 6; defensin beta 106A		

For research use only.

GENE DESCRIPTION	defensin, beta 106A
GENE ID	245909
ORIGIN SPECIES	Human
ACCESSIONS	BC100844
VECTOR	pDONR223
NCBI SYMBOL	DEFB106A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggacttctctcttctcttgcggtgctctcttctgacccagccaagaatgcatTTTTgatgagaatgcacaaactaaaggacatgcaagaacaattgcgggaaaaacgaagaacttattgctctctgccagaagtctgaaatgctgctggaccatccagccatgtgggagca
ttatagattaa