

DEFB106B-defensin, beta 106B Gene

Reference: PTXBC100846

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BD-6; DEFB-6; beta-defensin 106; beta-defensin 6; defensin, beta 106; defensin beta 106B		

For research use only.

GENE DESCRIPTION	defensin, beta 106B
GENE ID	503841
ORIGIN SPECIES	Human
ACCESSIONS	BC100846
VECTOR	pDONR223
NCBI SYMBOL	DEFB106B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggacttctctcttcttccgtgctctcttctgacccagccaagaatgcatTTTTgatgagaatgcacaaactaaaggacatgcaagaacaattgcgggaaaaatgaagaactattgctctctgccagaagtctctgaaatgctgctggaccatccagccatgtggagcat
tatagattaa