

DEFB104A-defensin, beta 104A Gene

Reference: PTXBC100848

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BD-4; DEFB-4; DEFB104; DEFB4; hBD-4; beta-defensin 104; defensin, beta 4; defensin beta 104A		

For research use only.

GENE DESCRIPTION	defensin, beta 104A
GENE ID	140596
ORIGIN SPECIES	Human
ACCESSIONS	BC100848
VECTOR	pDONR223
NCBI SYMBOL	DEFB104A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagagactgtgctgctattagccatttcttctactctatcaagatcttcagtgagaagcgaatttgaattggacagaatatgtggttatgggactgcccggtgccggaagaaatgtcgagccaagaatacagaattggaagatgtcccaacacctatgcctgtttgagaaaatggg
atgagagcttactgaatcgtacaaaacctga