

DEFB104B-defensin, beta 104B Gene

Reference: PTXBC100849

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BD-4; DEFB-4; hBD-4; beta-defensin 104; beta-defensin 4; defensin, beta 104; defensin beta 104B		

For research use only.

GENE DESCRIPTION	defensin, beta 104B
GENE ID	503618
ORIGIN SPECIES	Human
ACCESSIONS	BC100849
VECTOR	pDONR223
NCBI SYMBOL	DEFB104B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagagactgtgctgctattagccgtttctcttactctatcaagatctccagtgagaagcgaattggaattgacagaatatgtgttatggactgccgttgcggagaagaatgtcgcagccaagaatacagaattggaagatgtcccaacctatgcatgctgttgagaaatggatgagagcttactgaatcgtacaaaacctga