

IFNE1-interferon epsilon 1 Gene

Reference: PTXBC100873

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFNE1; IFN-E; IFNT1; INFE1; PRO655; interferon epsilon; IFN-epsilon; interferon epsilon 1; interferon tau-1		

For research use only.

GENE DESCRIPTION	interferon epsilon 1
GENE ID	338376
ORIGIN SPECIES	Human
ACCESSIONS	BC100873
VECTOR	pDONR223
NCBI SYMBOL	IFNE1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattatcaagcactcttgggaactgtgtggtgctgctgctcctaccactatctctctctagattgaaactgattatcttcagcaaaagacaagtgaatcaagaaagttaaactcttgaataagttgcaaaccttgcaattcagcagtgctaccacagggaaaaactttctctcctcag
aagtctttgagtcctcagcagtagcaaaaaaggacacactctggccattctccatgagatgcttcagcagatcttcagcctcttcagggaataattctctggtgtggaggaaaaccacacggagaaattcctcattcaacttcacacagctagaatacctagaagcactcatgggact
ggaagcagagaagctaagtgtggtggtggtgataaccttagattacaagttaaatgtacttccgaaggatccatgattacctggaaaaccaggactacagcactgtgcttggccattgtccaagtagaaatcagccgatgtctgttctgtgttcagtcacagaaaaactgagca
aacaaggagaccctgaacgacatgaagcaagagcttactacagagtttagaagcccgaggtag