

LYG2-lysozyme G-like 2 Gene

Reference: PTXBC100884

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LYGH; LYSG2; lysozyme g-like protein 2; lysozyme G-like 2; lysozyme g2		
For research use only.			

GENE DESCRIPTION	lysozyme G-like 2
GENE ID	254773
ORIGIN SPECIES	Human
ACCESSIONS	BC100884
VECTOR	pDONR223
NCBI SYMBOL	LYG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtatctccgtggtgttttgggactaattgcctcattggcacttcaggggtcctaccccttcagtcactcaatgaagcctcacctacatccacgctgtaccacggctgtatggggacatcatgacctgaagacctctggggcactgtgatgcaaacagtgtgatgaactgcgg
gatccgtggttctgaaatgtttgctgagatggatttgagggccataaaaccttaccagactctgatcaaaagatcgggagagacattgctggaccctgtctcatcgagccatcatctccagggaagccatggcggtatgtcctgcaagacgctgggaccacaggggacttaaat
ttggcttgatgcagcttgataaacaacgtaccaccctgtcgtgctgggtagcaaaagacaccttcacaggctactgggattctaacagagagaattaaggcaatccagaaaaattcccacgtggagtgttgcacgacaccaaagggtgtctctcagcttttaagtcaggaattga
agcgattgcccccacgacatagacaatgacttgcctcaatgatcattgctcagcctaagttctataaaagacaaagcttctag