



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD300E-CD300e molecule Gene

Reference: PTXBC100890

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD300e molecule; CD300e antigen; CD300LE; CLM-2; CLM2; CMRF35-A5; IREM-2; IREM2; PIgR-2; PIgR2; CMRF35-like molecule 2; CD300 antigen-like family member E; immune receptor expressed on myeloid cells 2; poly-Ig receptor 2; polymeric immunoglobulin receptor		

For research use only.

GENE DESCRIPTION	CD300e molecule
GENE ID	342510
ORIGIN SPECIES	Human
ACCESSIONS	BC100890
VECTOR	pDONR223
NCBI SYMBOL	CD300E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgctgctccagcctctactctctctcctcctctcaggctgtttgtctctgaaggcccccggctctgtgactggcactcgggggactctctgacagtgtgtgtctcgtatgagagcatgtacaaggatataacaagtactgtgtccgagagacagtacgacactcatgtgagacattgtggagacaaggagagaagagaaaggttggagagggaatggccgctgtccatcagagacaccggaggctctgccttctactgtgacatgtcagaacctcaatgaagatgtatgtctggatcttactgtgtcaaaatctagacagctgtgggtctctggattcatgttgcacgcgatccctcgacctgtatgggtgtatgtttcccccgaattacaaccccaaggaggaccacacatccaggccacacttccactctctctgtgtgtaacctctgcggaaacctcagaccagggaggtgtgtgaccaaaatcagggttccggctcagcagccctactctctctgtgtgtgtgtctctctgtgaagctg

More info: Please contact us for any question: contact@proteogenix.fr