

CD8B-CD8b molecule Gene

Reference: PTXBC100913

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD8b molecule; CD8B1; LEU2; LY3; LYT3; P37; T-cell surface glycoprotein CD8 beta chain; CD8 antigen, beta polypeptide 1 (p37); T lymphocyte surface glycoprotein beta chain		

For research use only.

GENE DESCRIPTION	CD8b molecule
GENE ID	926
ORIGIN SPECIES	Human
ACCESSIONS	BC100913
VECTOR	pDONR223
NCBI SYMBOL	CD8B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcggccgctgtggtctctctgcccgcagctgacagttctccatggcaactcagtcctccagcagaccctgcatacataaagggtgcaaaccaacaagatggtgatgctgtctctgcgaggctaaaatctccctcagtaacatgcgcactactggtgagacagcgcaggca
ccgagcagtgacagtcaccacgagttctctgcccctctgggattccgaaaagggaactatccacgggtgaagaggtggaacaggagaagatagctgtgttctgggatgcaagccggttcattctcaatctcacaagcgtgaagccggaagacagtgatctactctgcatgatctcggg
agccccgagctgaccttcgggaagggaactcagctgagtggtgtgattcttcccaccactgcccagcccaccaagaagtcacccctcaagaagagagtggtccggttaccagggccagagaccagaaggccacctttgtagccccatcaccttggcctgctgtggtggtggt
cctggttctgctggtttccctgggagtgccatccacctgtgctgcccggggagagagcccggttcgtttcatgaaacagcctcaagggaaggtatcatcaggaaacctttgtccccaatgcctgcatggatactacagcaataactacaacctcacagaagctctaaccatggatcct
gaaaacatag