

FGF14-fibroblast growth factor 14 Gene

Reference: PTXBC100920

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FGF-14; FHF-4; FHF4; SCA27; fibroblast growth factor 14; bA397O8.2; fibroblast growth factor homologous factor 4		

For research use only.

GENE DESCRIPTION	fibroblast growth factor 14
GENE ID	2259
ORIGIN SPECIES	Human
ACCESSIONS	BC100920
VECTOR	pDONR223
NCBI SYMBOL	FGF14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtaaaaccggtgccctcttcaggagaactgattcaaatattattatgcaaccacaaggatctcttcttcagggtgtctaaagctgctgattgcttttgcacaaatcaatgtggttctttggaacatttcagcaaaggaacgcatactgctcagtgcttttggcaagagtcttaagaa
aaacaagaacccaactgatccccagctcaagggtatagtgaccaggttatattgcaggcaaggctactacttgcaaatgcaccccgatggagctctcgatggaaccaaggatgacagcactaattctacactctcaacctcataccagtgaggactacgtgtgtgtccatccaggagtgaa
aacagggtgtatagccatgaatggagaaggttacctctacccatcagaacttttaccctgaatgcaagtttaagaatctgttttggaaaattattatgtaattactatccatgtgtacagacaacaggaatctgtagagcctggttttgggattaaataagggaaggcagctatgaaa
gggaacagagtaagaaacaaacaccagcagctattttaccacagccattggaagttgccatgtaccgagaaccatcttgcagatgattgggaaacggtcccgaagcctgggtgacgccagtaaaagcacaagtgctgtcgaataatgaatggaggcaaccagtcacaa
gagtaagacaacatag