

RASSF3-Ras association (RalGDS/AF-6) domain family member 3 Gene

Reference: PTXBC100952

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RASSF5; ras association domain-containing protein 3; Ras association (RalGDS/AF-6) domain family 3; Ras association (RalGDS/AF-6) domain family member 3; Ras association domain family member 3		

For research use only.

GENE DESCRIPTION	Ras association (RalGDS/AF-6) domain family member 3
GENE ID	283349
ORIGIN SPECIES	Human
ACCESSIONS	BC100952
VECTOR	pDONR223
NCBI SYMBOL	RASSF3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctgaattcaaatgggatttacactggcttcattaaagtacagatggaactctgcaaacctccacagacttctccaaattctggaaaactctctcccagtagcaatggctgtatgaatacacttcatacagcagcacaacactgtcgggaagtgtcggccctgtcaaaaagtctctgtgactgagagccctgccaaagtgtgacattataagcgtgtgcacagggaagacaaagtctacgctgcaagctctcagaccgggaacatccactctacctgcgtttggttagcagggccagaacagacacttagttttgtctctgtgaacatgaaattggagagtggaagccttcagcctccagaactacagaattcttgcgcatcttgacaaggaagaagatgaacagctgcagaacctgaagaggcgtacacagcctacaggcagaagctggaagaagccctccgtgaggtgtggaagcctgattaa