

C17orf65-chromosome 17 open reading frame 65 Gene

Reference: PTXBC100965

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf65; ASB16 antisense RNA 1		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 65
GENE ID	339201
ORIGIN SPECIES	Human
ACCESSIONS	BC100965
VECTOR	pDONR223
NCBI SYMBOL	C17orf65
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaggcctccagcgcacctcaggcgcactggctgcgcccgtagtcagcagtgccgcaaaaggacttcaggctccagttgggggagtcctggacggcctgcagcgcacagtcacatgggcgtgtggcccgcccatggggacctcagcgcggggcccgtaacgcag
cagcagctcggccagggcccccagccgttggcacaagcgttggcagcggcggtggtgcgttgcgccggccggccggccatcagctccagctccaggagccggcgccgcagcctggtgtgcctgcagctacctggccctcgcccccagcgcagccgtgttcagcgc
agtcctgcctgggtgtgcagggccacgtcggcccatgctccaggtacagagccacatgttgcctcaggccgcgcggccgcacgtgcagggcggtctcctggctctgcctgctccaggttcaccgtgcctcctgcttcagcagcaactggcgacctggggacaaga
gagcaaggtcctgaggaaaatcctgaaggcttgggacccctcagcttgcctgga