

LOC441251-Williams Beuren syndrome chromosome region 19 pseudogene Gene

Reference: PTXBC100974

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	Williams Beuren syndrome chromosome region 19 pseudogene
GENE ID	441251
ORIGIN SPECIES	Human
ACCESSIONS	BC100974
VECTOR	pDONR223
NCBI SYMBOL	LOC441251
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtggtggacaaatctgaggagtcgttggaggagccacggaaggtgtcgccttgagcctgaggagatctgggtggcggagatgctgtgtggcctcaagatgaagctgaagcgacggcgagtgctgctcctgagcaccagggcctcaacaggctgctt
gaggatcctgtcattaaaagattcctggcctgggacaaaggtctgagggtgtcggacaagtatctcctgctatggtcagatgtatttcagccggcgccctccctcctggcaataccaatgcattcatttctcctgctctctacctggccaatgacatgaggaggacgacgaggacc
ccaaacaaacatcttacttctgtatgggaagaccgctctcgcatacccttgcctcgttaagcgtcgggtccagttatgccgttgcatgaaccgagggccagggaagaccgctcagatagtcctgttccagaaacttcgggtccagttctctgttccatgagctgcaggcctgggttc
cccggaggagtgaggagatccaggccttatgacccagagcactgggtgtggcgcgagatcgctcgcttctctag