



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C12orf64-chromosome 12 open reading frame 64 Gene

Reference: PTXBC101016

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf64; DFNB84B; otogelin-like protein; otogelin like		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 64
GENE ID	283310
ORIGIN SPECIES	Human
ACCESSIONS	BC101016
VECTOR	pDONR223
NCBI SYMBOL	C12orf64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]