

## THSD3-thrombospondin, type I, domain containing 3 Gene

Reference: PTXBC101020

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | THSD3; isthmin-2; isthmin 2 homolog; thrombospondin and AMOP containing isthmin-like 1; thrombospondin and AMOP domain-containing isthmin-like protein 1; thrombospondin type-1 domain-containing protein 3; thrombospondin, type I, domain containing 3; isthm |                     |                |

For research use only.

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | thrombospondin, type I, domain containing 3                                                      |
| <b>GENE ID</b>          | 145501                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC101020                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | THSD3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgcgtgcgtccgcaccgagccgggctcctctctgctgctgctgctggcgccgctgctgagggcgccgctagggctccccgtgaagaagccgcggctccgcggaccacgcctgggagcctcacgaggtctcagcctccccagatcctagggcctctgaagga  
agaggaggaggcaccactgctccccagaaccacactgagggcagagccacacacatggatgctggactgctgagccagcagccatgacccaggcaacacacccctccaggacccagaggttactcgttgcggtgagctgcagaagctgcccggattggccaac  
acaacctgaglacccctaaccctgataccagggcttcagcctccccagatcctagggcctctgagggagaggaggagggcagcactgctccccagaaccacactgcaggcagagctacaccaacatggatgttgactgtcactgagccagcagcctgacccagggaatgccacg  
cctcccagaccagggaggttactccttctgctgagctgcagaagctgccagaattgttccacgcaacctttagtaccctaaacctgataaccaggtgacctcaaggtgtgaggagcccccaggccaggtgtcgtatagacctgttggctgagcccagcaatccccgccca  
ggatacccttagctggctgcccgcctctggtccttctctggggagactacaaggagagagaaaaagacagggccccagggggagagggggagggaaaaagaggaagacgaggactatccttcagaggatacaggggtgaggatcaagaggacaaagggaagatgaggaag  
agcaggcgctctggttcaatggaactacagacaactgggaccagggctggtggtccccggggattgggtctcaaggattctgctgactacgactatgagcctcagaaggagtgaggtccctggtctccctgagctgggaactgcagcactggcaagcagcagagggactgcccctgt  
ggctatggctgactgccaccagaccgtacctgtgactgcccctctgctgctggcactgaggacaaggacaccttggggcctccccagtgaggagtggaagctccttggcccgcaatgctacggacatgcatgaagatgtggacagctgtgaagatggctgaaactgcaagagcg  
acttctaatacagatctgagccagatgctgcgggacctgcccagctgcccgtgtgctaccactggagggcatggacagccctgtgagcctacaggagcagcaccagggccgagcttccgtggaggagatgccagtgccctcgcgagcgctggacatctaccagcccagg  
cgcgcttctgctgctgcttccatgctgctgctggggagagcagcacactggcccccagcactgctgctatgacgagagacagccggctgctgacccgtggcaagggcgccggcatgcccaacctcatcagcaccgacttctcacctaagctgacttcaagttcgacacgacccctggatc  
ctgtgcaagggggactggagcgcctccacgtgtgtcctcccaacaacggccgagcctgcaccgacaacccctggaggaggagtagctacacagttgcaggaggccaaggagtactag