



Reference: PTXBC101024

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	synaptotagmin-10; synaptotagmin X; sytX		
For research use only.			

GENE DESCRIPTION	synaptotagmin X
GENE ID	341359
ORIGIN SPECIES	Human
ACCESSIONS	BC101024
VECTOR	pDONR223
NCBI SYMBOL	SYT10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1