

## SASS6-spindle assembly 6 homolog (C. elegans) Gene

Reference: PTXBC101026

|                |                                                                                                     |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | MCPH14; SAS-6; SAS6; spindle assembly abnormal protein 6 homolog; SAS-6 centriolar assembly protein |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | spindle assembly 6 homolog (C. elegans)                                                          |
| <b>GENE ID</b>          | 163786                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC101026                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | SASS6                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgcaatcactagatgatctactaagcaactggactttacagaaagacattagcagaaaaaaacaagaattagataagttacggaatgaatggcgctcacatacagcagccttgacaacaagcattctcaggaaactgacaatgaaaaggaaaagccttcaggcacaggttcaat  
atcaacagcagcatgaacacagaaaaagatttagaaatcctccatcaacaaacatccaccagctacaaaacagactgtctgagttagaagtggttaataaagacttaaccgaaagaaataataaggagactccactattagagaacttaagcaaaacttctggtgtgagaggag  
ctacagcggactaagcaagaagtcctctcttgcgaagagagaattctactagatgttgaatgccagagaaagaaagcagcgttaatcagctacaacaaaagtgccagttttagaacaggaaatcaaggataaggaccagctgttttaagaacaaaaggcatttgatacaatcca  
ggaacaaaagggtgtttagaagaaaaggtgagaaaaatcaagtacaactaggaaagcttgaagctacataaaatcattatctgcagaactctgaaggcaaatgaaattatcaagaagttacaaggggatctgaaaactttaatgggttaagtgaattgaagaatacagttactattcagc  
aagaaaaactcttgctgagaaggagaaaaattacaaaaggaacaaaaggaaattacaagatgttggacagttcttcgaattaaagagcaagaggtatgcaattacaagaacaattagaagctacagttaaaaaacttgaagaaagcaaacacttctaaaaataatgaaaagttaac  
acgtggttaataaagaactaaatgaaatcagctagtgagaagcaagatgtattggaccttctactactccgctgcacattccagcagcaacacaatcagaagtggaatttctcctaactgaatgtggttgatggttagactgacttaccacactgtggattggttatcctgtctcctcg  
catttgattccagaatactctccctcattcgatctgccaataatcaccagccactgtgttcaggaaacaaaggttcagtttaatttgcagtttacaacccaatgcactactaggagatgttcagtcaggagcaactattatgctcctgtcactgataaggaaaatggtgaaaatgtagg  
gttggaaatccaaatcctgaagaaaagggaagatagcattcctttacgcggactcagccagaacctatttagtaattcagaccatcagagagatggcactttaggagcattacatactctccaacccacagcgctccctctgcgtcttcagcctatttccctgggcagttaccaaacagtta  
a