



Reference: PTXBC101048

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDH7L3; Cdh7; cadherin-20; cadherin 20, type 2; cadherin 20		
For research use only.			

GENE DESCRIPTION	cadherin 20, type 2
GENE ID	28316
ORIGIN SPECIES	Human
ACCESSIONS	BC101048
VECTOR	pDONR223
NCBI SYMBOL	CDH20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr