



Reference: PTXBC101063

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMYA4; CTRCT43; SMUNC45; UNC-45B; UNC45; protein unc-45 homolog B; cardiomyopathy associated 4; striated muscle UNC45; unc-45 homolog B; unc-45 myosin chaperone B		

For research use only.

GENE DESCRIPTION	unc-45 homolog B (C. elegans)
GENE ID	146862
ORIGIN SPECIES	Human
ACCESSIONS	BC101063
VECTOR	pDONR223
NCBI SYMBOL	UNC45B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr