



Reference: PTXBC101101

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAT2; TRAMD1; proton-coupled amino acid transporter 2; solute carrier family 36 (proton/amino acid symporter), member 2; tramdorin-1; solute carrier family 36 member 2		

For research use only.

GENE DESCRIPTION	solute carrier family 36 (proton/amino acid symporter), member 2
GENE ID	153201
ORIGIN SPECIES	Human
ACCESSIONS	BC101101
VECTOR	pDONR223
NCBI SYMBOL	SLC36A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]