

COX8C-cytochrome c oxidase subunit 8C Gene

Reference: PTXBC101126

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX8-3; cytochrome c oxidase subunit 8C, mitochondrial; COX VIII-3; cytochrome c oxidase polypeptide 8; cytochrome c oxidase polypeptide VIII; cytochrome c oxidase subunit 8-3; cytochrome c oxidase subunit VIII; cytochrome c oxidase subunit VIIIC; cytochr		

For research use only.

GENE DESCRIPTION	cytochrome c oxidase subunit 8C
GENE ID	341947
ORIGIN SPECIES	Human
ACCESSIONS	BC101126
VECTOR	pDONR223
NCBI SYMBOL	COX8C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgectctctgctggcgctgtctgccccgccactaccgcccgttggccctgctggcctgcagcccgctcccgcttgcacctcggggcccccgccagcgccctgtctgccgggaatggctgttggacttgtggtgtttttacgaccttctaacaccagctgcatat
gtgctaggcaacctgaagcagttcagaagggaattag