

QRFP-pyroglutamylated RFamide peptide Gene

Reference: PTXBC101127

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	prepro-QRFP; orexigenic neuropeptide QRFP; 26RFa; P518; P518 precursor protein; RF(Arg-Phe)amide family 26 amino acid peptide (P518); pyroglutamylated RFamide peptide		

For research use only.

GENE DESCRIPTION	pyroglutamylated RFamide peptide
GENE ID	347148
ORIGIN SPECIES	Human
ACCESSIONS	BC101127
VECTOR	pDONR223
NCBI SYMBOL	QRFP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtaaggccttaccctgatctacttctctctgcccgtggcgccctgcttcccttactggacagaagagagccacagacgcatgggtggcctcggagctggagaacgctggcgccgacctggccatggggcccgacccactccgtgtgggttcctctcgggtgctgaga
gcttcacagccacagccctgcttgcataccaggggtgagacatcggcgagagagcagctgctgctgcagattccgcttcgggagggagacgaaggcagtgaggccacggcttccctcctgctgcgggggagagaccagcgcccgtagggaaacctggctgaggagc
tcaatggctacagcaggaagaaggcggttcagcttcgcttcggtcggcggtga