

BAGE3-B melanoma antigen family, member 3 Gene

Reference: PTXBC101134

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT2.3; B melanoma antigen 3; cancer/testis antigen 2.3; cancer/testis antigen family 2, member 3; B melanoma antigen family member 3		

For research use only.

GENE DESCRIPTION	B melanoma antigen family, member 3
GENE ID	85318
ORIGIN SPECIES	Human
ACCESSIONS	BC101134
VECTOR	pDONR223
NCBI SYMBOL	BAGE3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggctggagtggtttttctgcatgtctgccagctgtccaagccaggtgatgaaggaggagtcctgtgtgagctggaggttgagcctgaagatggcacagctctcgatgtgcatgttgaagcactttggagccactatcaaatgctgtgaagagaaatgtaccagatgta
tcattatcctgtgtgctgaggagccgacagctttcaggaattcagtcacatctctgctttgtccagaacacattgaccaagctctgaaagatcgaaggaagatgcaaactgtgcagtgctgcacagccgggagacctcttag