

FERD3L-Fer3-like (Drosophila) Gene

Reference: PTXBC101135

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	N-TWIST; NATO3; NTWIST; PTFB; bHLHa31; fer3-like protein; Fer3-like; basic helix-loop-helix protein N-twist; class A basic helix-loop-helix protein 31; nephew of atonal 3; neuronal twist; pancreas-specific transcription factor b; Fer3 like bHLH transcript		

For research use only.

GENE DESCRIPTION	Fer3-like (Drosophila)
GENE ID	222894
ORIGIN SPECIES	Human
ACCESSIONS	BC101135
VECTOR	pDONR223
NCBI SYMBOL	FERD3L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccctatccggagagctcggtggacactacgggtgctggactcgtcgagacctgtccctggcctccccgagacgccctctctctgctgacttcgacccggggtctccttgggggaccagcccttgcgctccgagagggaagaccagaggatggcgcggtttgaagagg
gggacccagaagaagaggagtgcgaaagtggaccagggggacggagaagaggagggaagaagaggagcgcgggaagggtgtctccctattaggccgcccccaagaggaaggggtgatcacctacgccagcgccagggcccaacatccgcgaagggaagcggatgttca
acctcaacgaggcctttgaccagctgagggaagggtgccacgtttgcttacgagaaaagggtgtcccgatcgagacctccgctggccatgctctatatctcttcatgaccgagctcttgagagctgtgagaagaaggaaagcgctga